

PAW Manuka Wound Gel™



PAW Manuka Wound Gel™ is a sterile medical grade wound dressing made from *Leptospermum.sp* (Manuka) honey (80%) and natural oils and waxes for wound management. It is ideal for dogs, cats, horses and other companion animals with wounds, burns, grazes and as a general first aid.

PAW Manuka Wound Gel™ provides the following beneficial properties:

- Protects the wound
- Helps maintain natural wound pH
- Remove Malodour
- Reduces and prevents contamination, external bacteria colonisation and infection
- May assist in removal of slough and necrotic tissues
- Natural & non-toxic for animals
- Minimal heat treatment to preserve honey's natural properties
- Special blend of natural waxes and oils to reduce the stinging associated with pure honey
- Non-adherent gel to reduce trauma at dressing change.

Directions for Use

Wash around the wound with a sterile saline. Apply PAW Manuka Wound Gel™ liberally (approx. 3mm layer) directly to the wound or to a suitable dressing. In addition using a barrier cream can further protect the wound margin. Dressing change frequency will depend on clinical circumstances. It can vary from once per day for draining wounds to up to 7 days for unsoiled dressings. PAW Manuka Wound Gel™ can be removed from the wound by rinsing with sterile saline or water at dressing change.

Ingredients

The ingredients in PAW Manuka Wound Gel™ include Manuka Honey (80%), natural Gum, bees wax, and nourishing organic oils. This product is sterile and of medical grade.

Safety Directions

Not for use on heavily exudating wounds

Not recommended for body cavities

In the event of a slow or non-responding wound, please seek veterinary advice.

Dosage & Format

Available in a sterile 25g & 100g tube with induction seal and screw cap.

When to Use

Use on dogs, cats, horses and other companion animals with:

- Non-healing and chronic wounds
- Traumatic, acute and surgical wounds
- Malodorous and sloughy wounds
- Burns
- Abrasions & grazes
- Punctures, cuts & lacerations
- General first aid

Physical Description

A light brown gel formulation – resulting in reduced run off and easier dressing changes.